


Technical Datasheet

Rigid DL240 **Plant-Based**









Rigid DL240 Plant-Based is a high-performance rigid 3D printing resin which consists of 50% bio-based raw materials, offering a substantial reduction on net CO₂ emission compared to conventional resins. It is remarkably easy to handle and process, along with exhibiting outstanding properties.

Rigid DL240 Plant-Based has high accuracy, with over 98% of scanned data within +/- 100µm for dental models printed horizontally and over 83% of scanned data within +/- 100µm for dental models printed vertically, increasing output for overnight production. It enables quick design iterations by offering 250 and 350µm layer thickness print profiles.

Optimised for: • Fast & Accurate Prototyping

Dental Models for Aligner manufacturing

Unique features:








Exceptional surface finish and smooth feel



Quick design iterations



Rigid









Rigid DL240 Plant-Based Properties

| Tensile Properites | Green | Post-Cured | Method |
|------------------------------------|------------------------|------------------------------|------------|
| Tensile Modulus | 1210 MPa | 2440 MPa | ASTM D638 |
| Tensile Strength (Break) | 27 MPa | 56.6 MPa | ASTM D638 |
| Tensile Strength (Yield) | 24.3 MPa | 64.5 MPa | ASTM D638 |
| Elongation at Break | 19.2% | 6.1% | ASTM D638 |
| Flexural Properties | | | |
| Flexural Strength | - | 108 MPa | ASTM D790 |
| Flexural Modulus | - | 2656 MPa | ASTM D790 |
| Impact Properties | | | |
| Impact Strength Notched Izod | - | 12.2 J/m | ASTM D256 |
| Impact Strength Notched Izod | - | 2.2 kJ/m2 | ISO 180 |
| General Properties | | | |
| Shore Hardness | - | 88 Shore D | ASTM D2240 |
| HDT (@ 0.455 MPa) | - | 78.4°C | ASTM D648 |
| HDT (@ 1.82 MPa) | | 62.6°C | ASTM D648 |
| Water absorption (%)* after 24 hrs | - | 0.470% | ASTM D570 |
| Water absorption (%)* after 72 hrs | - | 0.625% | ASTM D570 |
| Water absorption (%)* after 7 days | - | 0.933% | ASTM D570 |
| Liquid Properties | Value | Method | |
| Viscosity | 580 cPs | At 25°C Brookfield spindle 3 | |
| Density | 1.10 g/cm ³ | - | |
| Storage | 10 <t>50°C</t> | _ | |

^{*} Post cured for 2 hours at 60°C with Photocentric Cure L2



Design & Print Orientation Consideration Parameters

There are some design guidelines for printing parts with Liquid Crystal Magna

| Properties | Parameters |
|------------------------------------|--|
| Minimum feature size (pins) | 0.5mm |
| Minimum hole diameter | 0.8mm |
| Minimum slot thickness | 0.4mm |
| Minimum wall thickness | 0.5mm |
| Overhangs | Successful for overhangs ≤15° |
| Round Dim Fit | Parts fit with no resistance at 0.01mm offset Click to view sample |
| Square Dim Fit | Parts fit perfectly with no resistance at 0.01mm offset Click to view sample |
| Minimum wall thickness unsupported | 0.5mm (add 0.25mm thicknes for every 10mm) |



Pre-Print Instructions

- 1. To print with Photocentric Liquid Crystal Magna, choose 'Rigid DL240Bio' at desired layer thickness when preparing your print file in Photocentric Studio.
- 2. Heat the resin to 30°C in the bottle.
- 3. Shake the resin bottle for 2 minutes before pouring into the resin vat.



Post-Print Instructions

- 1. Parts can be washed in 15 minutes using Photocentric Resin Cleaner or alternatively, in 10 minutes using Photocentric Resin Cleaner 30.
- 2. Make sure you do not exceed the recommended wash cycles as it might have adverse effect on the mechanical properties.
- 3. Once washed, rinse with warm water for 1-2 minutes
- 4. Dry with compressed air to remove any remaining water. Or alternatively, leave to air-dry.
- 5. To reach the ultimate mechanical properties: Place the platform into the Photocentric Cure L2 for a minimum of 2 hours at 60°C. for medium to large parts post curing might take between 3 to 4 hours.
 - If only 'dry to touch' finish is required, for example for dental models, 30 minutes post curing should be adequate.
- 6. Remove the platform from the Cure L2 and immediately leave it for 2 minutes under running cold water below 14°C for thermal shocking. Parts can be removed from the platform with minimal effort.





